

ARG70191 Mouse IGF2 recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Mouse IGF2 recombinant protein
Tested Application	SDS-PAGE
Target Name	IGF2
Species	Mouse
A.A. Sequence	Ala25 - Glu91
Expression System	E. coli
Activity	Active
Alternate Names	Somatomedin-A; Insulin-like growth factor II; IGF-II; PP9974; GRDF; C11orf43; T3M-11-derived growth factor

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images



ARG70191 Mouse IGF2 recombinant protein (Active) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70191 Mouse IGF2 recombinant protein (Active) (His-tagged, N-ter).