

ARG70193 Mouse M-CSF recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Mouse M-CSF recombinant protein
Tested Application	SDS-PAGE
Target Name	M-CSF
Species	Mouse
A.A. Sequence	Lys33 - Pro187
Expression System	E. coli
Activity	Active
Alternate Names	Macrophage colony-stimulating factor 1; CSF-1; Lanimostim; M-CSF; MCSF

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images



ARG70193 Mouse M-CSF recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70193 Mouse M-CSF recombinant protein (Active) (His-tagged, C-ter).