

ARG70200 Mouse beta NGF recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Mouse beta NGF recombinant protein
Tested Application	SDS-PAGE
Target Name	beta NGF
Species	Mouse
A.A. Sequence	Ser122 - Gly241
Expression System	E. coli
Activity	Active
Alternate Names	HSAN5; Beta-nerve growth factor; Beta-NGF; NGFB

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	20 mM sodium citrate and 0.2 M NaCl (pH 4.5)
Reconstitution	

Images



ARG70200 Mouse beta NGF recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70200 Mouse beta NGF recombinant protein (Active) (His-tagged, C-ter).