

ARG70202

Pig IL1 beta recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg

Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Pig IL1 beta recombinant protein
Tested Application	SDS-PAGE
Target Name	IL1 beta
Species	Pig
A.A. Sequence	Ala115 - Pro267
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce D10.G4.1 cells proliferation. The ED50 for this effect is < 3 ng/mL.
Alternate Names	Interleukin-1 beta; IL1-BETA; IL-1; IL-1 beta; Catabolin; IL1F2

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70202 Pig IL1 beta recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70202 Pig IL1 beta recombinant protein (Active) (His-tagged, C-ter).