

## ARG70203 Pig IL2 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Pig IL2 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL2
Species	Pig
A.A. Sequence	Ala21 - Thr154
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce proliferation in CTLL2 cells. The ED50 for this effect is < 0.5 ng/mL.
Alternate Names	TCGF; IL-2; lymphokine; Interleukin-2; Aldesleukin; T-cell growth factor

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

### Images



ARG70203 Pig IL2 recombinant protein (Active) (His-tagged, C-ter)  
SDS-PAGE image

SDS-PAGE analysis of ARG70203 Pig IL2 recombinant protein (Active) (His-tagged, C-ter).