

## ARG70204 Pig IL4 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Pig IL4 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL4
Species	Pig
A.A. Sequence	His25 - Cys133
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce proliferation in TF-1 cells. The ED50 for this effect is < 0.6 ng/mL.
Alternate Names	BSF-1; IL-4; B-cell stimulatory factor 1; BCGF1; Lymphocyte stimulatory factor 1; BCGF-1; Interleukin-4; Pitracinra; Binetrakin; BSF1

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

### Images



ARG70204 Pig IL4 recombinant protein (Active) (His-tagged, C-ter)  
SDS-PAGE image

SDS-PAGE analysis of ARG70204 Pig IL4 recombinant protein (Active)  
(His-tagged, C-ter).