

ARG70205 Pig IL6 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Pig IL6 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL6
Species	Pig
A.A. Sequence	Arg31 - Met212
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce proliferation in T1165.85.2.1 cells. The ED50 for this effect is < 1.5 ng/mL.
Alternate Names	B-cell stimulatory factor 2; CDF; HSF; BSF-2; Interferon beta-2; IL-6; IFNB2; CTL differentiation factor; Interleukin-6; HGF; Hybridoma growth factor; BSF2; IFN-beta-2

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70205 Pig IL6 recombinant protein (Active) (His-tagged, C-ter)
SDS-PAGE image

SDS-PAGE analysis of ARG70205 Pig IL6 recombinant protein (Active) (His-tagged, C-ter).