

ARG70207

Pig IL10 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg

Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Pig IL10 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL10
Species	Pig
A.A. Sequence	Ser19 - Asn175
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce proliferation in MC/9 β 2 cells. The ED50 for this effect is < 5 ng/mL.
Alternate Names	IL10A; GVHDS; TGIF; IL-10; Cytokine synthesis inhibitory factor; CSIF; Interleukin-10

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.01 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70207 Pig IL10 recombinant protein (Active) (His-tagged, C-ter)
SDS-PAGE image

SDS-PAGE analysis of ARG70207 Pig IL10 recombinant protein (Active) (His-tagged, C-ter).