

## ARG70208 Pig IL15 recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Pig IL15 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL15
Species	Pig
A.A. Sequence	Thr49 - Ser162
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce proliferation in TF-1 cells. The ED50 for this effect is < 5.5 ng/mL.
Alternate Names	IL-15; Interleukin-15

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.01 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

### Images



ARG70208 Pig IL15 recombinant protein (Active) (His-tagged, N-ter)  
SDS-PAGE image

SDS-PAGE analysis of ARG70208 Pig IL15 recombinant protein  
(Active) (His-tagged, N-ter).