

ARG70212 Pig GM-CSF recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Pig GM-CSF recombinant protein
Tested Application	SDS-PAGE
Target Name	GM-CSF
Species	Pig
A.A. Sequence	Ala18 - Lys144
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce proliferation in TF-1 cells. The ED50 for this effect is < 3 ng/mL.
Alternate Names	GM-CSF; Granulocyte-macrophage colony-stimulating factor; CSF; GMCSF; Molgramostin; Colony-stimulating factor; Sargramostim

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70212 Pig GM-CSF recombinant protein (Active) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70212 Pig GM-CSF recombinant protein (Active) (His-tagged, N-ter).