

ARG70216 SARS-CoV-2 nucleocapsid recombinant protein (His-SUMO tagged, N-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	HEK293 expressed, His-SUMO tagged (N-ter) SARS-CoV-2 nucleocapsid recombinant protein
Tested Application	SDS-PAGE
Target Name	SARS-CoV-2 nucleocapsid protein
Species	Virus
A.A. Sequence	Met1 - Ala419
Expression System	HEK293

Properties

Form	Powder
Purity	>98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



SARS-CoV-2 nucleocapsid recombinant protein (His-SUMO tagged, N-ter)

ARG70216 SARS-CoV-2 nucleocapsid recombinant protein (His-SUMO tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70216 SARS-CoV-2 nucleocapsid recombinant protein (His-SUMO tagged, N-ter).