

ARG70218 SARS-CoV-2 Spike recombinant protein (S1 Subunit) (His-tagged, C-ter)

Package: 100 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) SARS-CoV-2 Spike recombinant protein (S1 Subunit)
Tested Application	SDS-PAGE
Target Name	SARS-CoV-2 Spike protein (S1 Subunit)
Species	Virus
A.A. Sequence	Val16 - Cys671
Expression System	E. coli

Properties

Form	Powder
Purity	>98% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images



SARS-CoV-2 Spike
recombinant protein (S1 Subunit)

ARG70218 SARS-CoV-2 Spike recombinant protein (S1 Subunit) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70218 SARS-CoV-2 Spike recombinant protein (S1 Subunit) (His-tagged, C-ter).