

## ARG70226 Mouse IFN beta 1a recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

|                     |                                                                                    |
|---------------------|------------------------------------------------------------------------------------|
| Product Description | E. coli expressed, His-tagged (C-ter) Active Mouse IFN beta 1a recombinant protein |
| Tested Application  | SDS-PAGE                                                                           |
| Target Name         | IFN beta 1a                                                                        |
| Species             | Mouse                                                                              |
| A.A. Sequence       | Ile22 - Asn182                                                                     |
| Expression System   | E. coli                                                                            |
| Activity            | Active                                                                             |
| Alternate Names     | Fibroblast interferon; IFN-beta; IFNB; IFF; IFB; Interferon beta                   |

### Properties

|                   |                                                                               |
|-------------------|-------------------------------------------------------------------------------|
| Form              | Powder                                                                        |
| Purification Note | Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test. |
| Purity            | > 98% (by SDS-PAGE)                                                           |
| Buffer            | PBS (pH 7.4)                                                                  |
| Reconstitution    |                                                                               |

### Images



ARG70226 Mouse IFN beta 1a recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70226 Mouse IFN beta 1a recombinant protein (Active) (His-tagged, C-ter).