

ARG70284

Human TGFBR2 / TGF beta Receptor II recombinant protein (ECD) (Fc-His-tagged, C-ter)

Package: 50 µg

Store at: -20°C

Summary

Product Description	HEK293 expressed, Fc-His-tagged (C-ter) Human TGFBR2 / TGF beta Receptor II recombinant protein (ECD).
Tested Reactivity	Hu
Tested Application	Binding, SDS-PAGE
Target Name	TGFBR2 / TGF beta Receptor II (ECD)
Species	Human
A.A. Sequence	Thr23 - Asp159 of Human TGFBR2 / TGF beta Receptor II (NP_003233.4) with an Fc-6X His tag at the C-terminus.
Expression System	HEK293
Alternate Names	FAA3; AAT3; TbetaR-II; TGF-beta receptor type-2; LDS2B; MFS2; TGF-beta type II receptor; EC 2.7.11.30; TAAD2; TGFR-2; LDS2; Transforming growth factor-beta receptor type II; TGFbeta-RII; TGF-beta receptor type II; RIIC; LDS1B

Application Instructions

Application Note	Binding activity test: Measured by its binding ability in a functional ELISA. Immobilized recombinant human TGFBR1 at 2 µg/ml (100 µl/well) can bind recombinant human TGFBR2 with a linear range of 1-5 ng/ml.
------------------	---

Properties

Form	Powder
Purification Note	0.22 µm filter sterilized. Endotoxin level is 95% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	Reconstitute to a concentration of 0.1 - 0.5 mg/ml in sterile distilled water.
Storage instruction	For long term, lyophilized protein should be stored at -20°C or -80°C. After reconstitution, aliquot and store at -20°C for up to one month, at 2-8°C for up to one week. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	TGFBR2
Gene Full Name	transforming growth factor, beta receptor II (70/80kDa)
Function	Transmembrane serine/threonine kinase forming with the TGF-beta type I serine/threonine kinase receptor, TGFBR1, the non-promiscuous receptor for the TGF-beta cytokines TGFBR1, TGFBR2 and TGFBR3. Transduces the TGFBR1, TGFBR2 and TGFBR3 signal from the cell surface to the cytoplasm and is thus

regulating a plethora of physiological and pathological processes including cell cycle arrest in epithelial and hematopoietic cells, control of mesenchymal cell proliferation and differentiation, wound healing, extracellular matrix production, immunosuppression and carcinogenesis. The formation of the receptor complex composed of 2 TGFBR1 and 2 TGFBR2 molecules symmetrically bound to the cytokine dimer results in the phosphorylation and the activation of TGFBR1 by the constitutively active TGFBR2. Activated TGFBR1 phosphorylates SMAD2 which dissociates from the receptor and interacts with SMAD4. The SMAD2-SMAD4 complex is subsequently translocated to the nucleus where it modulates the transcription of the TGF-beta-regulated genes. This constitutes the canonical SMAD-dependent TGF-beta signaling cascade. Also involved in non-canonical, SMAD-independent TGF-beta signaling pathways. [UniProt]

Calculated Mw	65 kDa
PTM	Phosphorylated on a Ser/Thr residue in the cytoplasmic domain. [UniProt]
Cellular Localization	Cell membrane; Single-pass type I membrane protein. Membrane raft. [UniProt]

Images

