

ARG70358
Dengue virus NS1 recombinant protein (His-tagged, C-ter)Package: 5 µg
Store at: -20°C, -80°C

Summary

Product Description	HEK293 cell expressed, His-tagged (C-ter) Dengue virus NS1 recombinant protein
Tested Application	SDS-PAGE
Target Name	Dengue virus NS1 protein
Species	Virus
A.A. Sequence	Asp776 - Ala1127
Expression System	HEK293 cell
Alternate Names	Dengue virus NS1 protein; Dengue virus nonstructural glycoprotein NS1 protein

Properties

Form	Powder
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	