

## ARG70372 Mouse IFN gamma recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg  
Store at: -20°C

### Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Mouse IFN gamma recombinant protein
Tested Application	SDS-PAGE
Target Name	IFN gamma
Species	Mouse
A.A. Sequence	His23 - Cys155
Expression System	E. coli
Activity	Active
Alternate Names	IFNG; Interferon Gamma; Immune Interferon; IFN-Gamma; IMD; IFG; IFI

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

### Images



ARG70372 Mouse IFN gamma recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70372 Mouse IFN gamma recombinant protein (Active) (His-tagged, C-ter).