

ARG70377 Human IL17D recombinant protein (His-tagged, N-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (N-ter) Human IL17D recombinant protein
Tested Application	SDS-PAGE
Target Name	IL17D
Species	Human
A.A. Sequence	Ala18 - Pro202
Expression System	E. coli
Alternate Names	IL17D; Interleukin 17D; IL-17D; IL-27; Interleukin-17D; FLJ30846; IL27; Interleukin 27; Interleukin-27

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images



ARG70377 Human IL17D recombinant protein (His-tagged, N-ter)
SDS-PAGE image

SDS-PAGE analysis of ARG70377 Human IL17D recombinant protein (His-tagged, N-ter)