

**ARG70380**  
**Human Lymphotoxin Alpha / TNF beta recombinant protein (Active)**  
**(His-tagged, N-ter)**

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Human Lymphotoxin Alpha / TNF beta recombinant protein
Tested Application	SDS-PAGE
Target Name	Lymphotoxin alpha / TNF beta
Species	Human
A.A. Sequence	Leu1 - Leu171
Expression System	E. coli
Activity	Active
Alternate Names	LTA; Lymphotoxin Alpha; TNFSF1; TNFB; LT; Tumor Necrosis Factor Ligand Superfamily Member 1; Lymphotoxin-Alpha; LT-Alpha; TNF-Beta; Lymphotoxin Alpha (TNF Superfamily, Member 1); Tumor Necrosis Factor Ligand 1E; Tumor Necrosis Factor Beta; TNF Superfamily, Member 1; TNF Superfamily Member 1; TNLG1E

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	20 mM sodium citrate and 0.2 M NaCl (pH 3.5)
Reconstitution	

### Images



ARG70380 Human Lymphotoxin Alpha / TNF beta recombinant protein (Active) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70380 Human Lymphotoxin Alpha / TNF beta recombinant protein (Active) (His-tagged, N-ter)