

ARG70383 Human HMGB2 recombinant protein (His-SUMO tagged, N-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	HEK293 expressed, His-SUMO tagged (N-ter) Human HMGB2 recombinant protein
Tested Application	SDS-PAGE
Target Name	HMGB2
Species	Human
A.A. Sequence	Met1 - Pro187
Expression System	HEK293
Alternate Names	HMGB2; High Mobility Group Box 2; HMG2; High-Mobility Group (Nonhistone Chromosomal) Protein 2; High Mobility Group Protein B2; High Mobility Group Protein 2; HMG-2; High-Mobility Group Box 2

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70383 Human HMGB2 recombinant protein (His-SUMO tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70383 Human HMGB2 recombinant protein (His-SUMO tagged, N-ter)