

ARG70386 Human IFNW1 recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Human IFNW1 recombinant protein
Tested Application	SDS-PAGE
Target Name	IFNW1
Species	Human
A.A. Sequence	Cys24 - Ser195
Expression System	E. coli
Activity	Active
Alternate Names	IFNW1; Interferon Omega 1; IFN-Omega 1, Interferon Omega-1; Interferon Alpha-II-1; Interferon Omega-1

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70386 Human IFNW1 recombinant protein SDS-PAGE image

SDS-PAGE analysis of ARG70386 Human IFNW1 recombinant protein