

ARG70387 Human CD178 / Fas Ligand recombinant protein (Active) (His-SUMO tagged, N-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-SUMO tagged (N-ter) Active Human CD178 / Fas Ligand recombinant protein
Tested Application	SDS-PAGE
Target Name	CD178 / Fas Ligand
Species	Human
A.A. Sequence	Gln130 - Leu281
Expression System	E. coli
Activity	Active
Alternate Names	FASLG; Fas Ligand; APT1LG1; TNFSF6; CD178; FasL; Tumor Necrosis Factor Ligand Superfamily Member 6; Fas Ligand (TNF Superfamily, Member 6); Apoptosis Antigen Ligand; Fas Antigen Ligand; CD95 Ligand; CD95-L; CD95L; APTL; FASL; Tumor Necrosis Factor (Ligand) Superfamily, Member 6; Mutant Tumor Necrosis Factor Family Member 6; Apoptosis (APO-1) Antigen Ligand 1; Tumor Necrosis Factor Ligand 1A; CD178 Antigen; ALPS1B; TNLG1A

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images



ARG70387 Human CD178 / Fas Ligand recombinant protein (Active)
(His-SUMO tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70387 Human CD178 / Fas Ligand
recombinant protein (Active) (His-SUMO tagged, N-ter)