

ARG70390

Human GIF / Gastric intrinsic factor recombinant protein (His-tagged, C-ter)

Package: 100 µg, 20 µg

Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Human GIF / Gastric intrinsic factor recombinant protein
Tested Application	SDS-PAGE
Target Name	GIF / Gastric intrinsic factor
Species	Human
A.A. Sequence	Met1 - Tyr417
Expression System	E. coli
Alternate Names	CBLIF; Cobalamin Binding Intrinsic Factor; IFMH; INF; IF; TCN3; GIF; Gastric Intrinsic Factor (Vitamin B Synthesis); Gastric Intrinsic Factor; Intrinsic Factor

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images



ARG70390 Human GIF / Gastric intrinsic factor recombinant protein (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70390 Human GIF / Gastric intrinsic factor recombinant protein (His-tagged, C-ter)