

ARG70391 Human CHI3L1 recombinant protein (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Human CHI3L1 recombinant protein
Tested Application	SDS-PAGE
Target Name	CHI3L1
Species	Human
A.A. Sequence	Tyr22 - Thr383
Expression System	E. coli
Alternate Names	CHI3L1; Chitinase 3 Like 1; YKL40 YK-40; GP39; Chitinase 3-Like 1 (Cartilage Glycoprotein-39); Chitinase-3-Like Protein 1; Cartilage Glycoprotein 39; 39 KDa Synovial Protein; HCGP-39; CGP-39; YKL-40; GP-39; Cartilage Glycoprotein-39; HC-Gp39; HCGP-3P; YYL-40; ASRT7

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70391 Human CHI3L1 recombinant protein (His-tagged, C-ter)
SDS-PAGE image

SDS-PAGE analysis of ARG70391 Human CHI3L1 recombinant protein (His-tagged, C-ter)