

ARG70394 Mouse IL16 recombinant protein (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Mouse IL16 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL16
Species	Mouse
A.A. Sequence	His499 - Ser624
Expression System	E. coli
Alternate Names	IL16; Interleukin 16; PrIL-16; LCF; Lymphocyte Chemoattractant Factor; Pro-Interleukin-16; Prointerleukin 16; HsT19289; FLJ42735; FLJ16806; IL-16; Interleukin 16 (Lymphocyte Chemoattractant Factor); Neuronal Interleukin 16; PRIL16; NIL16

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70394 Mouse IL16 recombinant protein (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70394 Mouse IL16 recombinant protein (His-tagged, C-ter)