

ARG70395 Mouse IL17B recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Mouse IL17B recombinant protein
Tested Application	SDS-PAGE
Target Name	IL17B
Species	Mouse
A.A. Sequence	His21 - Phe180
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce IL-8 secretion in HepG2 cells. The ED50 for this effect is < 1.5 ng/mL.
Alternate Names	IL17B; Interleukin 17B; IL-17B; ZCYTO7; IL-20; NIRF; Neuronal Interleukin-17-Related Factor; Interleukin-17B; MGC138900; MGC138901; Neuronal Interleukin-17 Related Factor; Cytokine-Like Protein ZCYTO7; Interleukin-17 Beta; Cytokine Zcyto7; Interleukin 20; Interleukin-20; IL20

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70395 Mouse IL17B recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70395 Mouse IL17B recombinant protein (Active) (His-tagged, C-ter)