

ARG70396 Mouse IL17D recombinant protein (His-tagged, N-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (N-ter) Mouse IL17D recombinant protein
Tested Application	SDS-PAGE
Target Name	IL17D
Species	Mouse
A.A. Sequence	Ala25 - Arg205
Expression System	E. coli
Alternate Names	IL17D; Interleukin 17D; IL-17D; IL-27; Interleukin-17D; FLJ30846; IL27; Interleukin 27; Interleukin-27

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	20 mM sodium citrate and 0.2 M NaCl (pH 4.5)
Reconstitution	

Images



ARG70396 Mouse IL17D recombinant protein (His-tagged, N-ter)
SDS-PAGE image

SDS-PAGE analysis of ARG70396 Mouse IL17D recombinant protein (His-tagged, N-ter)