

ARG70397 Mouse IL31 recombinant protein (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Mouse IL31 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL31
Species	Mouse
A.A. Sequence	Thr24 - Cys163
Expression System	E. coli
Alternate Names	IL31; Interleukin 31; IL-31; Interleukin-31

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70397 Mouse IL31 recombinant protein (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70397 Mouse IL31 recombinant protein (His-tagged, C-ter)