

## ARG70398 Mouse IL38 recombinant protein (His-tagged, C-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (C-ter) Mouse IL38 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL38
Species	Mouse
A.A. Sequence	Met1 - Arg152
Expression System	E. coli
Alternate Names	IL1F10; Interleukin 1 Family Member 10; IL-1HY2; IL1-Theta; IL-1F10; FKSG75; Interleukin-1 Receptor Antagonist FKSG75; Interleukin-1 Family Member 10; Family Of Interleukin 1-Theta; Interleukin-1 HY2; Interleukin-38; FIL1 Theta; IL-1 Theta; MGC119832; MGC119833; MGC11983; IL1HY2; IL-38; Interleukin-1 Receptor Antagonist-Like FIL1 Theta; Interleukin 1 Family, Member 10 (Theta); Interleukin 1 Family Member 10 (Theta); IL-1F10 (Canonical Form IL-1F10a); Interleukin-1 Theta; FIL1- Theta; FIL1-Theta; FIL1T; IL38

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

### Images



ARG70398 Mouse IL38 recombinant protein (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70398 Mouse IL38 recombinant protein (His-tagged, C-ter)