

ARG70403 Mouse MIF recombinant protein (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Mouse MIF recombinant protein
Tested Application	SDS-PAGE
Target Name	MIF
Species	Mouse
A.A. Sequence	Met1 - Ala115
Expression System	E. coli
Alternate Names	MIF; Macrophage Migration Inhibitory Factor; GIF; Phenylpyruvate Tautomerase; GLIF; Macrophage Migration Inhibitory Factor (Glycosylation-Inhibiting Factor); Glycosylation-Inhibiting Factor; L-Dopachrome Tautomerase; L-Dopachrome Isomerase; MMIF; Epididymis Secretory Sperm Binding Protein; EC 5.3.3.12; EC 5.3.2.1

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70403 Mouse MIF recombinant protein (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70403 Mouse MIF recombinant protein (His-tagged, C-ter) recombinant protein