

ARG70406 Mouse CD255 / TWEAK recombinant protein (Active) (His-tagged, C-ter)


Package: 100 µg, 20 µg
Store at: -20°C

Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Mouse CD255 / TWEAK recombinant protein
Tested Application	SDS-PAGE
Target Name	CD255 / TWEAK
Species	Mouse
A.A. Sequence	His110 - His249
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce proliferation in HUVEC cells. The ED50 for this effect is < 0.2 µg/mL.
Alternate Names	TNFSF12; TNF Superfamily Member 12; TWEAK; DR3LG; APO3L; Tumor Necrosis Factor (Ligand) Superfamily, Member 12; Tumor Necrosis Factor Ligand Superfamily Member 12; APO3 Ligand; Tumor Necrosis Factor Superfamily Member 12; TNF-Related WEAK Inducer Of Apoptosis; TNF-Related Weak Inducer Of Apoptosis; Tumor Necrosis Factor Ligand 4A; APO3/DR3 Ligand; TNLG4A

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	



Mouse CD255 / TWEAK
recombinant protein

ARG70406 Mouse CD255 / TWEAK recombinant protein (Active)
(His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70406 Mouse CD255 / TWEAK
recombinant protein (Active) (His-tagged, C-ter)