

## ARG70408 Pig IL17 / IL17A recombinant protein (His-tagged, C-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (C-ter) Pig IL17 / IL17A recombinant protein
Tested Application	SDS-PAGE
Target Name	IL17 / IL17A
Species	Pig
A.A. Sequence	Gly26 - Ser155
Expression System	E. coli
Alternate Names	IL17A; Interleukin 17A; IL-17A; IL-17; CTLA8; IL17; Interleukin 17 (Cytotoxic T-Lymphocyte-Associated Serine Esterase 8); Cytotoxic T-Lymphocyte-Associated Protein 8; Cytotoxic T-Lymphocyte-Associated Antigen 8; Interleukin-17A; CTLA-8; Interleukin-17; ILA17

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

### Images



ARG70408 Pig IL17 / IL17A recombinant protein (His-tagged, C-ter)  
SDS-PAGE image

SDS-PAGE analysis of ARG70408 Pig IL17 / IL17A recombinant protein (His-tagged, C-ter)