

## ARG70409 Pig BMP4 recombinant protein (His-tagged, C-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

|                     |  |
|---------------------|--|
| Product Description | E. coli expressed, His-tagged (C-ter) Pig BMP4 recombinant protein   |
| Tested Application  | SDS-PAGE   |
| Target Name         | BMP4   |
| Species             | Pig  |
| A.A. Sequence       | Arg227 - Cys408  |
| Expression System   | E. coli  |
| Alternate Names     | BMP4; Bone Morphogenetic Protein 4; BMP2B; Bone Morphogenetic Protein 2B; BMP2B1; MCOPS6; BMP-2B; OFC11; BMP-4; ZYME; DVR4 |

### Properties

|                   |   |
|-------------------|---|
| Form              | Powder  |
| Purification Note | Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test. |
| Purity            | > 98% (by SDS-PAGE)   |
| Buffer            | 20 mM sodium citrate and 0.2 M NaCl (pH 3.5)                                  |
| Reconstitution    |   |

### Images



ARG70409 Pig Pig BMP4 recombinant protein (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of Pig BMP4 recombinant protein (His-tagged, C-ter)