

## ARG70410 Pig IFN gamma recombinant protein (Active) (His-tagged, C-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Pig IFNG recombinant protein
Tested Application	SDS-PAGE
Target Name	IFN gamma
Species	Pig
A.A. Sequence	Gln23 - Lys166
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to protect PK15 cells infected with encephalomyocarditis (EMC) virus. The ED50 for this effect is < 40 pg/mL.
Alternate Names	IFNG; Interferon Gamma; Immune Interferon; IFN-Gamma; IMD69; IFG; IFI

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.01 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

### Images



ARG70410 Pig IFN gamma recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70410 Pig IFN gamma recombinant protein (Active) (His-tagged, C-ter)