


ARG70411
Pig TNF alpha recombinant protein (Active) (His-tagged, C-ter)Package: 100 µg, 20 µg
Store at: -20°C

Summary

Product Description	E. coli expressed, His-tagged (C-ter) Active Pig TNF recombinant protein
Tested Application	SDS-PAGE
Target Name	TNF alpha
Species	Pig
A.A. Sequence	Leu77 - Leu232
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce cytotoxicity in PK15 cells in the presence of the actinomycin D. The ED50 for this effect is < 15 pg/mL.
Alternate Names	TNF; Tumor Necrosis Factor; TNF-Alpha; TNFSF2; TNFA; DIF; Tumor Necrosis Factor Ligand Superfamily Member 2; TNF-A; Tumor Necrosis Factor (TNF Superfamily, Member 2); Tumor Necrosis Factor Ligand 1F; Tumor Necrosis Factor-Alpha; Tumor Necrotic Factor Alph; TNF Superfamily, Member 2; TNF, Macrophage-Derived; TNF, Monocyte-Derived; APC1 Protein Cachectin; TNLG1F

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	



ARG70411 Pig TNF alpha recombinant protein (Active) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70411 Pig TNF alpha recombinant protein (Active) (His-tagged, C-ter)