

ARG70416 Mouse HDGF recombinant protein (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Mouse HDGF recombinant protein
Tested Application	SDS-PAGE
Target Name	HDGF
Species	Mouse
A.A. Sequence	Met1 - Leu237
Expression System	E. coli
Alternate Names	HMG1L2; HMG-1L2; HDGF; Heparin Binding Growth Factor; HMG1L2; High Mobility Group Protein 1-Like 2; Hepatoma Derived Growth Factor

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70416 Mouse HDGF recombinant protein (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70416 Mouse HDGF recombinant protein (His-tagged, C-ter)