

ARG70417 Pig VEGF recombinant protein (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Pig VEGFA recombinant protein
Tested Application	SDS-PAGE
Target Name	VEGF
Species	Pig
A.A. Sequence	Ala27 - Arg190
Expression System	E. coli
Alternate Names	VEGFA; Vascular Endothelial Growth Factor A; VPF; VEGF-A; VEGF; Vascular Permeability Factor

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images



ARG70417 Pig VEGF recombinant protein (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70417 Pig VEGF recombinant protein (His-tagged, C-ter)