

ARG70418 Pig IGF1 recombinant protein (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Pig IGF1 recombinant protein
Tested Application	SDS-PAGE
Target Name	IGF1
Species	Pig
A.A. Sequence	Gly49 - Ala118
Expression System	E. coli
Alternate Names	IGF1; Insulin Like Growth Factor 1; IGF-I; IGF; IGF; Insulin-Like Growth Factor 1 (Somatomedin C); Insulin-Like Growth Factor I ; Mechano Growth Factor; Somatomedin-C; IGF1A; MGF; Insulin-Like Growth Factor IB; Somatomedin C; IBP1

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images



ARG70418 Pig IGF1 recombinant protein (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70418 Pig IGF1 recombinant protein (His-tagged, C-ter)