

ARG70419 Pig EGF recombinant protein (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Pig EGF recombinant protein
Tested Application	SDS-PAGE
Target Name	EGF
Species	Pig
A.A. Sequence	Asn970 - Arg1022
Expression System	E. coli
Alternate Names	EGF; Epidermal Growth Factor; Pro-Epidermal Growth Factor; Epidermal Growth Factor (Beta-Urogastrone); Beta-Urogastrone; HOMG4; URG

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images



ARG70419 Pig EGF recombinant protein (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70419 Pig EGF recombinant protein (His-tagged, C-ter)