

## ARG70420 Human CCL2 / MCP1 recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Human CCL2 / MCP1 recombinant protein
Tested Application	SDS-PAGE
Target Name	CCL2 / MCP1
Species	Human
A.A. Sequence	Gln24 - Thr99
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to chemoattract BaF3 cells transfected with human CCR2A. The ED50 for this effect is < 20 ng/mL.
Alternate Names	CCL2; C-C Motif Chemokine Ligand 2; MCP-1; MCP1; MCAF; HC11; Monocyte Chemotactic And Activating Factor Monocyte Secretory Protein JE; SMC-CF; GDCF-2; SCYA2

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

### Images



ARG70420 Human CCL2 / MCP1 recombinant protein (Active) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of Human CCL2 / MCP1 recombinant protein (Active) (His-tagged, N-ter)