

ARG70433

Mouse CCL3 / MIP1 alpha recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg


Store at: -20°C

Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Mouse CCL3 / MIP1 alpha recombinant protein
Tested Application	SDS-PAGE
Target Name	CCL3 / MIP1 alpha
Species	Mouse
A.A. Sequence	Ala24 - Ala92
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to chemoattract human PBMCs using a concentration range of 10.0 - 100.0 ng/mL. Note: Results may vary from different PBMC donors.
Alternate Names	CCL3; C-C Motif Chemokine Ligand 3; MIP-1-Alpha; G0S19-1; LD78ALPHA; SCYA3; Small Inducible Cytokine A3 (Homologous To Mouse Mip-1a); Macrophage Inflammatory Protein 1-Alpha

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	



ARG70433 Mouse CCL3 / MIP1 alpha recombinant protein (Active)
(His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70433 Mouse CCL3 / MIP1 alpha
recombinant protein (Active) (His-tagged, N-ter)