

ARG70434

Mouse CCL4 / MIP1 beta recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg


Store at: -20°C

Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Mouse CCL4 / MIP1 beta recombinant protein
Tested Application	SDS-PAGE
Target Name	CCL4 / MIP1 beta
Species	Mouse
A.A. Sequence	Ala24 - Asn92
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to chemoattract human PBMCs using a concentration range of 20.0 - 200.0 ng/mL. Note: Results may vary from different PBMC donors.
Alternate Names	CCL4; C-C Motif Chemokine Ligand 4; MIP-1-Beta; AT744.1; SCYA4; LAG1; Small Inducible Cytokine A4 (Homologous To Mouse Mip-1b); Macrophage Inflammatory Protein 1-Beta; Lymphocyte Activation Gene 1 Protein; G-26 T-Lymphocyte-Secreted Protein; Chemokine (C-C Motif) Ligand 4; T-Cell Activation Protein 2; C-C Motif Chemokine 4; MIP-1-Beta(1-69); SIS-Gamma

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	



ARG70434 Mouse CCL4 / MIP1 beta recombinant protein (Active)
(His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70434 Mouse CCL4 / MIP1 beta
recombinant protein (Active) (His-tagged, N-ter)