

ARG70436

Mouse CXCL2 / MIP2 recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg

Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Mouse CXCL2 / MIP2 recombinant protein
Tested Application	SDS-PAGE
Target Name	CXCL2 / MIP2
Species	Mouse
A.A. Sequence	Ala28 - Asn110
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to chemoattract BaF3 cells transfected with human CXCR2. The ED50 for this effect is < 0.5 ng/mL.
Alternate Names	CXCL2; C-X-C Motif Chemokine Ligand 2; SCYB2; CINC-2a; MIP-2a; MGSA-B; GROb; GRO2; Macrophage Inflammatory Protein 2-Alpha; Chemokine (C-X-C Motif) Ligand 2

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70436 Mouse CXCL2 / MIP2 recombinant protein (Active) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70436 Mouse CXCL2 / MIP2 recombinant protein (Active) (His-tagged, N-ter)