

ARG70437

Mouse CXCL3 / GRO gamma recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg


Store at: -20°C

Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Mouse CXCL3 recombinant protein
Tested Application	SDS-PAGE
Target Name	CXCL3 / GRO gamma
Species	Mouse
A.A. Sequence	Ala28 - Gln87
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to chemoattract BaF3 cells transfected with human CXCR2. The ED50 for this effect is < 80 ng/mL.
Alternate Names	CXCL3; C-X-C Motif Chemokine Ligand 3; SCYB3; CINC-2b; MIP-2b; GROg; GRO3; Macrophage Inflammatory Protein 2-Beta; Chemokine (C-X-C Motif) Ligand 3; Growth-Regulated Protein Gamma

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	



ARG70437 Mouse CXCL3 / GRO gamma recombinant protein (Active)
(His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70437 Mouse CXCL3 / GRO gamma
recombinant protein (Active) (His-tagged, N-ter)