


ARG70439
Mouse CXCL7 / NAP2 (40-113) recombinant protein (Active) (His-tagged, N-ter)Package: 100 µg, 20 µg
Store at: -20°C

Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Mouse CXCL7 / NAP2 recombinant protein
Tested Application	SDS-PAGE
Target Name	CXCL7 / NAP2
Species	Mouse
A.A. Sequence	Lys40 - Tyr113
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to chemoattract BaF3 cells transfected with human CXCR2. The ED50 for this effect is < 1 µg/mL.
Alternate Names	PPBP; Pro-Platelet Basic Protein; SCYB7; CTAP3; CXCL7; MDGF; LDGF; TGB1; PBP; Platelet Basic Protein; Beta-TG; CTAPIII; LA-PF4; THBGB1; B-TG1; NAP-2; TGB; Connective Tissue-Activating Peptide III; Chemokine (C-X-C Motif) Ligand 7

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	



ARG70439 Mouse CXCL7 / NAP2 (40-113) recombinant protein
(Active) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70439 Mouse CXCL7 / NAP2 (40-113)
recombinant protein (Active) (His-tagged, N-ter)