

## ARG70444 Mouse CXCL16 recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Mouse CXCL16 recombinant protein
Tested Application	SDS-PAGE
Target Name	CXCL16
Species	Mouse
A.A. Sequence	Asn27 - Pro114
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to chemoattract BaF3 cells transfected with mouse CXCR6. The ED50 for this effect is < 3 ng/mL.
Alternate Names	CXCL16; C X C Motif Chemokine Ligand 16; SR PSOX; SRPSOX; CXCLG16; Scavenger Receptor For Phosphatidylserine And Oxidized Low Density Lipoprotein

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

### Images



ARG70444 Mouse CXCL16 recombinant protein (Active) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70444 Mouse CXCL16 recombinant protein (Active) (His-tagged, N-ter)