

ARG70446 Pig CXCL9 / MIG recombinant protein (His-tagged, N-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (N-ter) Pig CXCL9 / MIG recombinant protein
Tested Application	SDS-PAGE
Target Name	CXCL9 / MIG
Species	Pig
A.A. Sequence	Thr23 - Thr126
Expression System	E. coli
Alternate Names	CXCL9; C-X-C Motif Chemokine Ligand 9; SCYB9; Crg-10; Humig; CMK; MIG

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70446 Pig CXCL9 / MIG recombinant protein (His-tagged, N-ter)
SDS-PAGE image

SDS-PAGE analysis of ARG70446 Pig CXCL9 / MIG recombinant protein (His-tagged, N-ter).