

ARG70447

Pig CXCL10 / IP10 recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg

Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Pig CXCL10 / IP10 recombinant protein
Tested Application	SDS-PAGE
Target Name	CXCL10 / IP10
Species	Pig
A.A. Sequence	Pro23 - Ala104
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to chemoattract BaF3 cells transfected with human CXCR2. The ED50 for this effect is < 4 ng/mL.
Alternate Names	CXCL10; C-X-C Motif Chemokine Ligand 10; IP-10; GIP-10; SCYB10; IFI10 ; Crg-2; Mob-1; INP10; C7

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70447 Pig CXCL10 / IP10 recombinant protein (Active) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70447 Pig CXCL10 / IP10 recombinant protein (Active) (His-tagged, N-ter)