

## ARG70448 Pig CXCL11 / I-TAC recombinant protein (His-tagged, N-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (N-ter) Pig CXCL11 / I-TAC recombinant protein
Tested Application	SDS-PAGE
Target Name	CXCL11 / I-TAC
Species	Pig
A.A. Sequence	Phe22 - Val100
Expression System	E. coli
Alternate Names	CXCL11; C-X-C Motif Chemokine Ligand 11; I-TAC; H174; IP-9; SCYB11; SCYB9B; B-R1; Small Inducible Cytokine Subfamily B (Cys-X-Cys), Member 11; Interferon-Inducible T-Cell Alpha Chemoattractant

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

### Images



ARG70448 Pig CXCL11 / I-TAC recombinant protein (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70448 Pig CXCL11 / I-TAC recombinant protein (His-tagged, N-ter)