

ARG70449 Pig CXCL3 / GRO gamma recombinant protein (His-tagged, N-ter)

Package: 100 µg, 20 µg
Store at: -20°C

Summary

Product Description	E. coli expressed, His-tagged (N-ter) Pig CXCL3 recombinant protein
Tested Application	SDS-PAGE
Target Name	CXCL3 / GRO gamma
Species	Pig
A.A. Sequence	Val24 - Ala111
Expression System	E. coli
Alternate Names	CXCL3; C-X-C Motif Chemokine Ligand 3; SCYB3; CINC-2b; MIP-2b; GROg; GRO3; Macrophage Inflammatory Protein 2-Beta; Chemokine (C-X-C Motif) Ligand 3; Growth-Regulated Protein Gamma

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	