

## ARG70450 Human MMP2 recombinant protein (His-tagged, C-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (C-ter) Human MMP2 recombinant protein
Tested Application	SDS-PAGE
Target Name	MMP2
Species	Human
A.A. Sequence	Tyr120 - Cys660
Expression System	E. coli
Alternate Names	MMP2; Matrix Metalloproteinase 2; TBE-1; CLG4A; 72 kDa Type IV Collagenase; Matrix Metalloproteinase-2; EC 3.4.24.24; MMP-2; CLG4; Matrix Metalloproteinase 2 (Gelatinase A, 72kDa Gelatinase, 72kDa Type IV Collagenase)

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

### Images



ARG70450 Human MMP2 recombinant protein (His-tagged, C-ter)  
SDS-PAGE image

SDS-PAGE analysis of ARG70450 Human MMP2 recombinant protein  
(His-tagged, C-ter)